MEDIC-2



BIOSYSTEM

Body Infrared Thermometer Accuracy: +/- 0.2°C



INFRARED FOREHEAD THERMOMETER 32.0 °C-43.0 °C



Forehead temperature measurement mode

0 °C-100 °C



Object temperature measurement mode



## INFRARED FOREHEAD THERMOMETER MEDIC-2



TEMPERATURE MEASUREMENT MODE

Tap the settings button to switch

INFRARED FOREHEAD THERMOMETER

MEDIC-2





**BEEP ALERT** 

to indicate high temperature



### INFRARED FOREHEAD THERMOMETER MEDIC-2



#### **BACKLIGHT**

For visibility in the dark

### NON-CONTACT

For optimal hygiene

1-5CM

range



32

Memory Storage



**2**s

Turn on



**30**s

Automatic Shut-down







# INFRARED FOREHEAD THERMOMETER MEDIC-2



#### Features of product usage

Accurate	The measurement error is < 0.3℃.
Fast	The measuring time is 0.5 second.
Ease-to-use	Instant measurement, easy to operate.
Non-contact	Measure forehead, no need to contact skin, safe and healthy, therefore no cross infection.
Screen show	Large-screen LCD shows clear reading even in dark night.
Multiple-use	It can measure human body temperature, after switching to object surface temperature model with measurement range of 0-100°C, it can measure room temperature, and temperature of baby bottle, nce paste and bath water.
Auto-shutdown	Auto-shutdown in 7 seconds, saving energy.
High temperature alarm	Prompt high temperature by buzzer, and the buzzer can be set to shutdown.